

For Research Use Only

CEACAM8 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68683-1-PBS



Basic Information

Catalog Number:

68683-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC026263

GeneID (NCBI):

1088

UNIPROT ID:

P31997

Full Name:

carcinoembryonic antigen-related
cell adhesion molecule 8

Calculated MW:

38 kDa

Purification Method:

Protein G purification

CloneNo.:

2B4F1

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50039-1, CEACAM8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68683-1-PBS. Detection antibody: 68683-2-PBS. Standard: Eg32133. Range: 46.88-3000 pg/mL.