

For Research Use Only

TNFSF18 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68684-2-PBS



Basic Information

Catalog Number:

68684-2-PBS

Concentration:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC069319

GeneID (NCBI):

8995

UNIPROT ID:

Q9UNG2

Full Name:

tumor necrosis factor (ligand)
superfamily, member 18

Calculated MW:

177 aa, 20 kDa

Purification Method:

Protein A purification

CloneNo.:

1C7F5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50040-1, TNFSF18 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68684-1-PBS. Detection antibody: 68684-2-PBS. Standard: Eg32094. Range: 62.5-4000 pg/mL.