

For Research Use Only

# AFP Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68686-2-PBS



## Basic Information

Catalog Number:

68686-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC027881

GeneID (NCBI):

174

UNIPROT ID:

P02771

Full Name:

alpha-fetoprotein

Calculated MW:

69 kDa

Purification Method:

Protein G purification

CloneNo.:

1D11G5

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50072-1, AFP  
Monoclonal Matched Antibody Pair, PBS Only.  
Capture antibody: 68686-1-PBS. Detection  
antibody: 68686-2-PBS. Standard: Eg0539. Range:  
31.25-2000 pg/mL