

For Research Use Only

S100A14 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68690-2-PBS



Basic Information

Catalog Number: 68690-2-PBS	GenBank Accession Number: BC005019	Purification Method: Protein A purification
Concentration: 1mg/ml	GeneID (NCBI): 57402	CloneNo.: 1F9H11
Source: Mouse	UNIPROT ID: Q9HCY8	
Isotype: IgG1	Full Name: S100 calcium binding protein A14	
Immunogen Catalog Number: AG29970	Calculated MW: 12 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50082-1, S100A14 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68690-1-PBS. Detection antibody: 68690-2-PBS. Standard: Ag29970. Range: 0.156-10 ng/mL



Serum of sixteen individual healthy human donors was measured. The human S100A14 concentration of detected samples was determined to be 4.20 ng/mL with a range of 0.31 - 24.71 ng/mL