

For Research Use Only

Transferrin/TF Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68697-2-PBS



Basic Information

Catalog Number:

68697-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

EG0114

GenBank Accession Number:

BC059367

GeneID (NCBI):

7018

UNIPROT ID:

P02787

Full Name:

transferrin

Calculated MW:

698 aa, 77 kDa

Purification Method:

Protein A purification

CloneNo.:

1G5E1

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Sandwich ELISA standard curve of MP50076-1, Transferrin/TF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68697-1-PBS. Detection antibody: 68697-2-PBS. Standard: Eg0114. Range: 23.4-1500 pg/mL.



Urine of eight individual healthy human donors was measured. The Transferrin/TF concentration of detected samples was determined to be 597.9 ng/mL with a range of 102.5-2,011.6 ng/mL.



Serum of sixteen individual healthy human donors was measured. The Transferrin/TF concentration of detected samples was determined to be 1.9 mg/mL with a range of 1.0-4.5 mg/mL.