

For Research Use Only

# CEACAM6 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68699-1-PBS



## Basic Information

Catalog Number:

68699-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC005008

GeneID (NCBI):

4680

UNIPROT ID:

P40199

Full Name:

carcinoembryonic antigen-related  
cell adhesion molecule 6 (non-  
specific cross reacting antigen)

Calculated MW:

37 kDa

Purification Method:

Protein G purification

CloneNo.:

2B4E1

## Applications

Tested Applications:

IF/ICC, Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Sandwich ELISA standard curve of MP50071-1, CEACAM6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68699-1-PBS. Detection antibody: 68699-2-PBS. Standard: Eg0056. Range: 31.25-2000 pg/mL