

For Research Use Only

CEACAM6 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68699-2-PBS



Basic Information

Catalog Number:

68699-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

BC005008

GeneID (NCBI):

4680

UNIPROT ID:

P40199

Full Name:

carcinoembryonic antigen-related
cell adhesion molecule 6 (non-
specific cross reacting antigen)

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

2D7G5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50071-1, CEACAM6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68699-1-PBS. Detection antibody: 68699-2-PBS. Standard: Eg0056. Range: 31.25-2000 pg/mL