

For Research Use Only

CD200 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68700-1-PBS



Basic Information

Catalog Number:

68700-1-PBS

Concentration:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

EG0006

GenBank Accession Number:

BC022522

GeneID (NCBI):

4345

UNIPROT ID:

P41217

Full Name:

CD200 molecule

Calculated MW:

269 aa, 30 kDa

Purification Method:

Protein G purification

CloneNo.:

1B1C10

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50062-1, CD200 Monoclonal Matched Antibody Pair - PBS only. 68700-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0006. 68700-2-PBS was HRP conjugated as the detection antibody. Range: 31.3-2000 pg/mL.



Serum of four individual healthy human donors was measured. The human CD200 concentration of detected samples was determined to be 374.2 pg/mL with a range of 164.2 - 490.7 pg/mL.