

For Research Use Only

CD200 Monoclonal antibody, PBS Only

Catalog Number: 68700-3-PBS



Basic Information

Catalog Number: 68700-3-PBS	GenBank Accession Number: BC022522	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 4345	CloneNo.: 2G8D1
Source: Mouse	UNIPROT ID: P41217	
Isotype: IgG1	Full Name: CD200 molecule	
Immunogen Catalog Number: EG0006	Calculated MW: 269 aa, 30 kDa	

Applications

Tested Applications:
Indirect ELISA
Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data