

For Research Use Only

GM-CSF Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68702-1-PBS



Basic Information

Catalog Number:

68702-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

HZ-1002

GenBank Accession Number:

GeneID (NCBI):

1437

Full Name:

colony stimulating factor 2
(granulocyte-macrophage)

Purification Method:

Protein A purification

CloneNo.:

2A5B1

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50067-1, GM-CSF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 69003-1-PBS. Detection antibody: 68702-1-PBS. Standard: Eg0189. Range: 7.8-500 pg/mL.