

For Research Use Only

Leptin R Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68703-2-PBS



Basic Information

Catalog Number:

68703-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG29621

GenBank Accession Number:

NM_002303

GeneID (NCBI):

3953

UNIPROT ID:

P48357

Full Name:

leptin receptor

Calculated MW:

132 kDa

Purification Method:

Protein A purification

CloneNo.:

2G6G6

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50066-1, LEPR
Monoclonal Matched Antibody Pair, PBS Only.
Capture antibody: 68703-1-PBS. Detection
antibody: 68703-2-PBS. Standard: Ag29621. Range:
39.06-2500 pg/mL