

For Research Use Only

CD46 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68707-1-PBS



Basic Information

Catalog Number: 68707-1-PBS	GenBank Accession Number: BC030594	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 4179	CloneNo.: 1A8A1
Source: Mouse	UNIPROT ID: P15529	
Isotype: IgG1	Full Name: CD46 molecule, complement regulatory protein	
Immunogen Catalog Number: AG28581	Calculated MW: 384 aa, 43 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50041-1, CD46 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68707-1-PBS. Detection antibody: 68707-2-PBS. Standard: Ag28581. Range: 31.25-1000 pg/mL