

For Research Use Only

CD46 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68707-2-PBS



Basic Information

Catalog Number:

68707-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28581

GenBank Accession Number:

BC030594

GeneID (NCBI):

4179

UNIPROT ID:

P15529

Full Name:

CD46 molecule, complement regulatory protein

Calculated MW:

384 aa, 43 kDa

Purification Method:

Protein G purification

CloneNo.:

2C6D12

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50041-1, CD46 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68707-1-PBS. Detection antibody: 68707-2-PBS. Standard: Ag28581. Range: 31.25-1000 pg/mL