

For Research Use Only

NASP Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68710-2-PBS



Basic Information

Catalog Number: 68710-2-PBS	GenBank Accession Number: BC010105	Purification Method: Protein G purification
Source: Mouse	GeneID (NCBI): 4678	CloneNo.: 1F6H2
Isotype: IgG1	UNIPROT ID: P49321	
Immunogen Catalog Number: AG1824	Full Name: nuclear autoantigenic sperm protein (histone-binding)	
	Calculated MW: 788 aa, 85 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51848-1, NASP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68710-1-PBS. Detection antibody: 68710-2-PBS.



Cytometric bead array standard curve of MP51848-1, NASP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68710-1-PBS. Detection antibody: 68710-2-PBS. Standard: Ag1824. Range: 0.781-100 ng/mL.