

For Research Use Only

CEP135 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68715-2-PBS



Basic Information

Catalog Number:

68715-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG19122

GenBank Accession Number:

BC136535

GeneID (NCBI):

9662

UNIPROT ID:

Q66GS9

Full Name:

centrosomal protein 135kDa

Calculated MW:

1140 aa, 134 kDa

Purification Method:

Protein G purification

CloneNo.:

2G9F2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51850-1, CEP135 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68715-1-PBS. Detection antibody: 68715-2-PBS.



Cytometric bead array standard curve of MP51850-1, CEP135 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68715-1-PBS. Detection antibody: 68715-2-PBS. Standard:Ag19122. Range: 0.781-100 ng/mL