

For Research Use Only

SPDEF Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68716-2-PBS



Basic Information

Catalog Number:

68716-2-PBS

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG33739

GenBank Accession Number:

BC021299

GeneID (NCBI):

25803

UNIPROT ID:

O95238

Full Name:

SAM pointed domain containing ets transcription factor

Calculated MW:

335 aa, 38 kDa

Purification Method:

Protein G purification

CloneNo.:

2B4E9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51889-1, SPDEF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68716-2-PBS. Detection antibody: 68716-1-PBS.



Cytometric bead array standard curve of MP51889-1, SPDEF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68716-2-PBS. Detection antibody: 68716-1-PBS. Standard:Ag33739. Range: 0.781-100 ng/mL.