

For Research Use Only

KIAA0319L Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68717-2-PBS



Basic Information

Catalog Number: 68717-2-PBS	GenBank Accession Number: BC031672	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 79932	CloneNo.: 2D10C6
Source: Mouse	UNIPROT ID: Q8I2AO	
Isotype: IgG1	Full Name: KIAA0319-like	
Immunogen Catalog Number: AG15263	Calculated MW: 1049 aa, 116 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50926-1, KIAA0319L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68717-2-PBS. Detection antibody: 68717-3-PBS. Standard: Ag15263. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50926-3, KIAA0319L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68717-2-PBS. Detection antibody: 68717-6-PBS. Standard: Ag15263. Range: 0.781-100 ng/mL