

For Research Use Only

KIAA0319L Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68717-3-PBS



Basic Information

Catalog Number:

68717-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG15263

GenBank Accession Number:

BC031672

GeneID (NCBI):

79932

UNIPROT ID:

Q8IZAO

Full Name:

KIAA0319-like

Calculated MW:

1049 aa, 116 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2B10G12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50926-1, KIAA0319L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68717-2-PBS. Detection antibody: 68717-3-PBS. Standard:Ag15263. Range: 0.098-100 ng/mL