

For Research Use Only

ALPL Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68719-3-PBS



Basic Information

Catalog Number:

68719-3-PBS

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG33678

GenBank Accession Number:

BC021289

GeneID (NCBI):

249

UNIPROT ID:

P05186

Full Name:

alkaline phosphatase,
liver/bone/kidney

Calculated MW:

57 kDa

Purification Method:

Protein A purification

CloneNo.:

2C10H7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array sample test of MP51482-1, ALPL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68719-2-PBS. Detection antibody: 68719-3-PBS.



Cytometric bead array standard curve of MP51482-1, ALPL Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68719-2-PBS. Detection antibody: 68719-3-PBS. Standard: Ag33678. Range: 0.781-100 ng/mL.