

For Research Use Only

SAMHD1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68732-1-PBS

Basic Information

Catalog Number: 68732-1-PBS	GenBank Accession Number: BC036450	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 25939	CloneNo.: 1D10H7
Source: Mouse	UNIPROT ID: Q9Y3Z3	
Isotype: IgG1	Full Name: SAM domain and HD domain 1	
Immunogen Catalog Number: AG18766	Calculated MW: 626 aa, 72 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve ofMP50020-1, SAMHD1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68732-1-PBS. Detection antibody: HRP-conjugated 68732-2-PBS. Standard: Ag18766. Range: 0.313-20 ng/mL.