

For Research Use Only

# TFF1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68733-2-PBS

## Basic Information

Catalog Number:

68733-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4541

GenBank Accession Number:

BC032811

GeneID (NCBI):

7031

UNIPROT ID:

P04155

Full Name:

trefoil factor 1

Calculated MW:

84 aa, 9 kDa

Purification Method:

Protein G purification

CloneNo.:

2A5D3

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Sandwich ELISA standard curve of MP50027-1, TFF1  
Monoclonal Matched Antibody Pair, PBS Only.  
Capture antibody: 68733-1-PBS. Detection  
antibody: HRP-conjugated 68733-2-PBS. Standard:  
Ag4541. Range: 0.625-20 ng/mL