

For Research Use Only

ATP1A1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68735-1-PBS

Basic Information

Catalog Number:

68735-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5990

GenBank Accession Number:

BC050359

GeneID (NCBI):

476

ENSEMBL Gene ID:

ENSG00000163399

UNIPROT ID:

P05023

Full Name:

ATPase, Na⁺/K⁺ transporting, alpha 1 polypeptide

Calculated MW:

113 kDa

Purification Method:

Protein G purification

CloneNo.:

2A2A11

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50088-1, ATP1A1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68735-1-PBS. Detection antibody: HRP-conjugated 68735-2-PBS. Standard: Ag5990. Range: 0.195-50 ng/mL.