

For Research Use Only

PD-L1/CD274 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68745-1-PBS



Basic Information

Catalog Number:

68745-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC074984

GeneID (NCBI):

29126

UNIPROT ID:

Q9NZQ7

Full Name:

CD274 molecule

Calculated MW:

290 aa, 33 kDa

Purification Method:

Protein G purification

CloneNo.:

2C12F12

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50090-1, PD-L1/CD274 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68745-1-PBS. Detection antibody: HRP-conjugated 68745-2-PBS. Standard:Eg31414. Range:0.049-1.563 ng/mL.