

For Research Use Only

# CD80 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68746-1-PBS



## Basic Information

Catalog Number:

68746-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM\_005191

GeneID (NCBI):

941

UNIPROT ID:

P33681

Full Name:

CD80 molecule

Purification Method:

Protein G purification

CloneNo.:

1A10C6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50098-1, CD80/B7-1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68746-1-PBS. Detection antibody: 68746-2-PBS. Standard:Eg0060. Range: 0.391-100 ng/mL