

For Research Use Only

# FGFR4 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68749-1-PBS

## Basic Information

Catalog Number:

68749-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC011847

GeneID (NCBI):

2264

UNIPROT ID:

P22455

Full Name:

fibroblast growth factor receptor 4

Calculated MW:

88 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3B3D7

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50104-1, FGFR4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68749-1-PBS. Detection antibody: 68749-2-PBS. Standard:Eg0221. Range: 1.563-100 ng/mL.