

For Research Use Only

# FIGF Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68750-2-PBS



## Basic Information

Catalog Number:

68750-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC027948

GeneID (NCBI):

2277

UNIPROT ID:

O43915

Full Name:

c-fos induced growth factor (vascular endothelial growth factor D)

Calculated MW:

40 kDa

Purification Method:

Protein G purification

CloneNo.:

1B1F4

## Applications

Tested Applications:

Cytometric bead array, Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50105-2, FIGF Monoclonal Matched Antibody Pair, PBS Only. 68750-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0184. 68750-2-PBS was HRP conjugated as the detection antibody. Range: 62.5-4000 pg/mL



Cytometric bead array standard curve of MP50105-1, FIGF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68750-1-PBS. Detection antibody: 68750-2-PBS. Standard:Eg0184. Range: 3.125-100 ng/mL