

For Research Use Only

FIGF Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68750-3-PBS

Basic Information

Catalog Number:

68750-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC027948

GeneID (NCBI):

2277

UNIPROT ID:

O43915

Full Name:

c-fos induced growth factor (vascular
endothelial growth factor D)

Calculated MW:

40 kDa

Purification Method:

Protein G purification

CloneNo.:

1E10C11

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50105-2, FIGF Monoclonal Matched Antibody Pair, PBS Only. 68750-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0184. 68750-2-PBS was HRP conjugated as the detection antibody. Range: 62.5-4000 pg/mL.