

For Research Use Only

# CTLA-4/CD152 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68751-3-PBS



## Basic Information

Catalog Number:

68751-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC070162

GeneID (NCBI):

1493

UNIPROT ID:

P16410

Full Name:

cytotoxic T-lymphocyte-associated  
protein 4

Calculated MW:

223 aa, 25 kDa

Purification Method:

Protein G purification

CloneNo.:

1F5E7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50106-1, CTLA4/CD152 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68751-3-PBS. Detection antibody: 68751-4-PBS. Standard:Eg0018. Range: 3.125-100 ng/mL