

For Research Use Only

TIGIT Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68752-1-PBS

Basic Information

Catalog Number:

68752-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM_173799.4

GeneID (NCBI):

201633

UNIPROT ID:

Q495A1-1

Full Name:

T cell immunoreceptor with Ig and ITIM domains

Calculated MW:

26 kDa

Purification Method:

Protein G purification

CloneNo.:

2D2C7

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50107-1, TIGIT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68752-1-PBS. Detection antibody: 68752-2-PBS. Standard:Eg0748. Range: 3.125-100 ng/mL.