

For Research Use Only

# IL-1 alpha Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68757-1-PBS



## Basic Information

Catalog Number:

68757-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC013142

GeneID (NCBI):

3552

UNIPROT ID:

P01583

Full Name:

interleukin 1, alpha

Calculated MW:

271 aa, 31 kDa

Purification Method:

Protein G purification

CloneNo.:

3A1G3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50113-1, IL1A/IL-1 alpha Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68757-1-PBS. Detection antibody: 68758-1-PBS. Standard:HZ-1320. Range: 0.098-3.125 ng/mL