

For Research Use Only

VEGFR1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68762-2-PBS

Basic Information

Catalog Number:

68762-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

NM_002019

GeneID (NCBI):

2321

UNIPROT ID:

P17948-1

Full Name:

fms-related tyrosine kinase 1
(vascular endothelial growth
factor/vascular permeability factor
receptor)

Calculated MW:

151kd

Purification Method:

Protein G purification

CloneNo.:

2B4C8

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50118-1, VEGFR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68762-1-PBS. Detection antibody: 68762-2-PBS. Standard:Eg0438. Range: 6.25-100 ng/mL