

For Research Use Only

# TMEM16A/DOG1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68764-1-PBS



## Basic Information

Catalog Number:

68764-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG3320

GenBank Accession Number:

BC027590

GeneID (NCBI):

55107

UNIPROT ID:

Q5XXA6

Full Name:

anoctamin 1, calcium activated  
chloride channel

Calculated MW:

986 aa, 114 kDa

Purification Method:

Protein A purification

CloneNo.:

1D11G7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50123-1, ANO1/TMEM16A/DOG1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68764-1-PBS. Detection antibody: 68764-2-PBS. Standard:Ag3320. Range: 0.098-100 ng/mL.