

For Research Use Only

# EIF6 Monoclonal antibody, PBS Only

Catalog Number: 68765-3-PBS



## Basic Information

Catalog Number:

68765-3-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG30249

GenBank Accession Number:

BC001119

GeneID (NCBI):

3692

UNIPROT ID:

P56537

Full Name:

eukaryotic translation initiation factor 6

Calculated MW:

27 kDa

Purification Method:

Protein A purification

CloneNo.:

1F2A8

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68765-3-Ig (EIF6 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68765-3-PBS in a different storage buffer formulation.