

For Research Use Only

ANKRD13B Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68766-2-PBS

Basic Information

Catalog Number:

68766-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG22725

GenBank Accession Number:

BC032554

GeneID (NCBI):

124930

UNIPROT ID:

Q86YJ7

Full Name:

ankyrin repeat domain 13B

Calculated MW:

626 aa, 70 kDa

Purification Method:

Protein G purification

CloneNo.:

1A4B3

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50125-1, ANKRD13B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68766-1-PBS. Detection antibody: 68766-2-PBS. Standard:Ag22725. Range: 1.563-100 ng/mL