

For Research Use Only

IL-1 beta Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68767-2-PBS



Basic Information

Catalog Number:

68767-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC008678

GeneID (NCBI):

3553

ENSEMBL Gene ID:

ENSG00000125538

UNIPROT ID:

P01584

Full Name:

interleukin 1, beta

Calculated MW:

269 aa, 31 kDa

Purification Method:

Protein G purification

CloneNo.:

3H8A3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP80001-3, IL-1 beta Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 82696-3-PBS. Detection antibody: 68767-2-PBS. Standard:HZ-1164. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP80001-4, IL-1 beta Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 82696-5-RR. Detection antibody: 68767-2-PBS. Standard:HZ-1164. Range:3.125-100 ng/mL.