

For Research Use Only

GDF-15 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68768-5-PBS



Basic Information

Catalog Number: 68768-5-PBS	GenBank Accession Number: BC008962	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 9518	CloneNo.: 2H6B3
Source: Mouse	UNIPROT ID: Q99988	
Isotype: IgG1	Full Name: growth differentiation factor 15	
Immunogen Catalog Number: AG33750	Calculated MW: 34 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50133-8, GDF15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68768-5-PBS. Detection antibody: HRP-conjugated 68768-1-PBS. Standard:Ag33750. Range: 1.563-100 ng/mL.



Cytometric bead array of LnCap cells supernatant using MP50133-8, GDF-15 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68768-5-PBS. Detection antibody: 68768-1-PBS. Complete culture medium was served as blank control.