

For Research Use Only

Pancreatic Lipase Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68770-1-PBS



Basic Information

Catalog Number:

68770-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27928

GenBank Accession Number:

BC014309

GeneID (NCBI):

5406

UNIPROT ID:

P16233

Full Name:

pancreatic lipase

Calculated MW:

52 kDa

Purification Method:

Protein G purification

CloneNo.:

2A7A8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50137-1, PNLIP/Pancreatic Lipase Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68770-1-PBS. Detection antibody: 68770-2-PBS. Standard:Ag27928. Range: 3.125-100 ng/mL.