

For Research Use Only

# TNFR2 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68771-1-PBS

## Basic Information

Catalog Number:

68771-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30497

GenBank Accession Number:

BC052977

GeneID (NCBI):

7133

UNIPROT ID:

P20333

Full Name:

tumor necrosis factor receptor  
superfamily, member 1B

Calculated MW:

48 kDa

Purification Method:

Protein G purification

CloneNo.:

3B5F4

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP50139-1, TNFR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68771-1-PBS. Detection antibody: HRP-conjugated 68771-2-PBS. Standard:Ag30497. Range: 0.098-25 ng/mL.