

For Research Use Only

TNFR2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68771-2-PBS

Basic Information

Catalog Number:

68771-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG30497

GenBank Accession Number:

BC052977

GeneID (NCBI):

7133

UNIPROT ID:

P20333

Full Name:

tumor necrosis factor receptor
superfamily, member 1B

Calculated MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

3D6E8

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50139-1, TNFR2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68771-1-PBS. Detection antibody: HRP-conjugated 68771-2-PBS. Standard: Ag30497. Range: 0.098-25 ng/mL.