

For Research Use Only

BMP-7 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68772-2-PBS



Basic Information

Catalog Number:

68772-2-PBS

Concentration:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

HZ-1229

GenBank Accession Number:

GeneID (NCBI):

655

Full Name:

bone morphogenetic protein 7

Purification Method:

Protein G purification

CloneNo.:

3E10E2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array of cell line supernatant using MP50138-1, BMP-7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68772-1-PBS. Detection antibody: 68772-2-PBS. Supernatant (50 µ L/well) from the following cell lines were assayed: MCF-7, LNCaP.



Cytometric bead array standard curve of MP50138-1, BMP-7 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68772-1-PBS. Detection antibody: 68772-2-PBS. Standard:HZ-1229. Range: 3.125-100 ng/mL.