

For Research Use Only

CALCA/Calcitonin Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68774-3-PBS



Basic Information

Catalog Number:

68774-3-PBS

Concentration:

1000 μ g/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG26141

GenBank Accession Number:

BC093753

GeneID (NCBI):

796

UNIPROT ID:

P01258

Full Name:

Calcitonin

Calculated MW:

141 aa, 15 kDa

Purification Method:

Protein G purification

CloneNo.:

1B8C4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array in cell lysate using MP50293-2, CALCA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68891-1-PBS. Detection antibody: 68774-3-PBS. Cell lysate: TT, HepG2.



Cytometric bead array standard curve of MP50293-2, CALCA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68891-1-PBS. Detection antibody: 68774-3-PBS. Standard: Ag26141. Range: 0.781-100 ng/mL.