

For Research Use Only

PPP2R1A Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68778-1-PBS



Basic Information

Catalog Number: 68778-1-PBS	GenBank Accession Number: BC001537	Purification Method: Protein G purification
Concentration: 1mg/ml	GeneID (NCBI): 5518	CloneNo.: 1D1C1
Source: Mouse	UNIPROT ID: P30153	
Isotype: IgG1	Full Name: protein phosphatase 2 (formerly 2A), regulatory subunit A, alpha isoform	
Immunogen Catalog Number: AG8661	Calculated MW: 589 aa, 65 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50144-1, PPP2R1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68778-1-PBS. Detection antibody: 68778-2-PBS. Standard:Ag8661. Range: 6.25-100 ng/mL