

For Research Use Only

PGDS Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68783-4-PBS



Basic Information

Catalog Number:

68783-4-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG34483

GenBank Accession Number:

BC005939

GeneID (NCBI):

5730

UNIPROT ID:

P41222

Full Name:

prostaglandin D2 synthase 21kDa (brain)

Calculated MW:

21 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2B3E1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sample test of MP50148-2, PGDS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68783-4-PBS. Detection antibody: 68783-5-PBS.



Cytometric bead array standard curve of MP50148-2, PGDS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68783-4-PBS. Detection antibody: 68783-5-PBS. Standard:Ag34483. Range: 0.781-100 ng/mL.