

For Research Use Only

ACTR8 Monoclonal antibody

Catalog Number: 68784-1-Ig



Basic Information

Catalog Number:

68784-1-Ig

Concentration:

1000 ug/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG34088

GenBank Accession Number:

BC032744

GeneID (NCBI):

93973

UNIPROT ID:

Q9H981

Full Name:

ARP8 actin-related protein 8 homolog (yeast)

Calculated MW:

329 aa, 37 kDa

Observed MW:

70 kDa

Purification Method:

Protein G purification

CloneNo.:

2B3B2

Recommended Dilutions:

WB 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293T cells, LNCaP cells, Jurkat cells, K-562 cells, Daudi cells, HSC-T6 cells, NIH/3T3 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68784-1-Ig (ACTR8 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



LNCaP cells were subjected to SDS PAGE followed by western blot with 68784-1-Ig (ACTR8 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.