

For Research Use Only

# SLC7A11 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68785-1-PBS



## Basic Information

Catalog Number:

68785-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30611

GenBank Accession Number:

BC012087

GeneID (NCBI):

23657

UNIPROT ID:

Q9UPY5

Full Name:

solute carrier family 7, (cationic amino acid transporter, y+ system) member 11

Calculated MW:

55 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3D7H10

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50151-1, SLC7A11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68785-1-PBS. Detection antibody: 68785-2-PBS. Standard:Ag30611. Range: 1.563-100 ng/mL.