

For Research Use Only

RIPK3 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68786-1-PBS

Basic Information

Catalog Number:

68786-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG11759

GenBank Accession Number:

BC062584

GeneID (NCBI):

11035

UNIPROT ID:

Q9Y572

Full Name:

receptor-interacting serine-threonine
kinase 3

Calculated MW:

518 aa, 57 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A10E4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50152-1, RIPK3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68786-1-PBS. Detection antibody: 68786-2-PBS. Standard:Ag11759. Range: 0.781-100 ng/mL.