

For Research Use Only

# Cyclin D1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68790-1-PBS

## Basic Information

Catalog Number:

68790-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25384

GenBank Accession Number:

BC000076

GeneID (NCBI):

595

UNIPROT ID:

P24385

Full Name:

cyclin D1

Calculated MW:

295 aa, 34 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G1H6

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Sandwich ELISA standard curve of MP80003-1, Cyclin D1 Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 68790-1-PBS. Detection antibody: 82681-1-PBS. Standard: Ag25384. Range: 0.188-12 ng/mL