

For Research Use Only

Galectin-3 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68791-1-PBS

Basic Information

Catalog Number:

68791-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC001120

GeneID (NCBI):

3958

UNIPROT ID:

P17931

Full Name:

lectin, galactoside-binding, soluble, 3

Calculated MW:

26 kDa

Purification Method:

Protein G purification

CloneNo.:

1G1B2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50157-1, Galectin 3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60207-3-PBS. Detection antibody: 68791-1-PBS. Standard:Eg0227. Range: 1.563-100 ng/mL