

For Research Use Only

# NME1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68792-2-PBS



## Basic Information

Catalog Number: 68792-2-PBS	GenBank Accession Number: BC018994	Purification Method: Protein G purification
Concentration: 1 mg/ml	GeneID (NCBI): 4830	CloneNo.: 2C6C7
Source: Mouse	UNIPROT ID: P15531	
Isotype: IgG1	Full Name: non-metastatic cells 1, protein (NM23A) expressed in	
Immunogen Catalog Number: AG28529	Calculated MW: 17 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP50159-1, NME1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68792-1-PBS. Detection antibody: 68792-2-PBS. Standard:Ag28529. Range: 0.391-100 ng/mL.