

For Research Use Only

# CEBP Alpha/CEBPA Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68793-1-PBS



## Basic Information

Catalog Number:

68793-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG34094

GenBank Accession Number:

NM\_004364

GeneID (NCBI):

1050

UNIPROT ID:

P49715

Full Name:

CCAAT/enhancer binding protein (C/EBP), alpha

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

3C10G4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP50160-1, CEBP Alpha/CEBPA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68793-1-PBS. Detection antibody: 68793-2-PBS. Standard:Ag34094. Range: 1.563-50 ng/mL.